

Saliva Collection Kits

OMNI offers two non-invasive, easy to use options for collecting saliva samples for SARS-CoV-2 testing.

OMNI International's Saliva Collection Kits provide an easy to use, non-invasive alternative for sample collection when testing for SARS-CoV-2. These sample collection kits are used to collect, store and transport saliva DNA samples and are ideal for use in subsequent procedures such as PCR.

Saliva Collection Kit with Cup

Saliva sample is discharged into the collection cup, then transferred using the disposable pipette into the 2 mL collection tube.

Kit Contents:

- Bio Hazard Bag
- 2 mL Collection Tube
- Collection Cup
- Disposable Pipette
- Alcohol Wipe




Saliva Collection Kit with Funnel

Saliva sample is discharged directly into the 2 mL collection tube via the provided funnel.

Kit Contents:

- Bio Hazard Bag
- 2 mL Collection Tube
- Funnel
- Alcohol Wipe



Ordering Information

Kit	Quantity	Order #	Price
Saliva Collection Kit with Funnel	100 per pack	26-030	\$475 USD
Saliva Collection Kit with Cup	100 per pack	26-031	\$300 USD

For research use only. Not for use in diagnostic procedures.



935 Cobb Place Blvd.
Kennesaw, GA 30144
800.776.4431 • 770.421.0058
Sales@omni-inc.com
www.omni-inc.com

Rev.A